# **Genetic test for** identification of parents (parentage testing)

### Lecture



Aleš Knoll, Petra Bartoňová

Mendel University in Brno, Faculty of AgriScience





### **Basic terms**

parentage: parenthood

paternity: fatherhood

Individual identification: determining that the DNA sample belongs to a specific individual (molecular fingerprinting)

Genetic type: a set of genotypes (alleles) characteristic and unmistakable for everyone, also known as *genetic profile* 



### Parentage testing and individual identification in animals

- importance verification of origin and identification
- livestock in CZ (it varies by country)
  - pig: determining of the genetic type (identification) of boars in breeding farms
  - cattle: genetic type before selection for breeding (bull breeding), verification and certification of the origin of breeding bulls
  - horse: verification and certification of the origin of foals after insemination (always for English thoroughbred, even if there is no insemination)
  - other species: dog, exotic birds, protected animals CITES (origin and identification)



### **Importance**

Importance in animal genetics, forensic genetics and medicine

#### Animals:

frequent errors in records incorrect parentage data to estimate breeding value you buy an animal that does not belong to the declared (quality) parents

#### Human:

decisive

paternity disputes identification of the offender, or victims transplantation and other medical procedures where kinship is

ISAGREED



### Development and use of methods

Previously: blood group polymorphism and protein

New: DNA variability (DNA fingerprinting)

- VNTR (variable number of tandem repeats, minisatellites)

#### Currently:

- STR (short tandem repeats, microsatellites)
- SNP (latest, increasingly used)

**VNTR** 

#### Requirements:

- high reliability of rejection and (or) confirmation of parentage
- relatively easy and quick laboratory testing
- high reproducibility
- the lowest price





### Microsatellites

- short tandem repetitions composed of mono, di, tri or tetra nucleotide motifs
- highly polymorphous ⇒ different alleles (5-20)

e.g. 
$$(GC)_n$$
 n = 5-25 repetitions

- lenght polymorphism i.e. if the MS sequence is amplified by PCR, we detect by electrophoresis a different fragment size (allele) due to a different number of repetitions
- the allele designation is therefore the base pair size



### **Mikrosatelites**

(GC)<sub>n</sub>

n=10

156 bp

↓ primer forward mikrosatelite primer reverse ↓

5 CATTGAATCGGTATCAT(N)<sub>60</sub>AGCGCGCGCGCGCGCGCT(N)<sub>40</sub>ACGTTAATTCATTCGTG

152 bp

5 CATTGAATCGGTATCAT(N)<sub>60</sub>AGCGCGCT(N)<sub>40</sub>ACGTTAATTCATTCGTG


142 bp

<sup>5</sup>CATTGAATCGGTATCAT(N)<sub>60</sub>AGCGCGCGCGCGCGCGCGCGCT(N)<sub>40</sub>ACGTTAATTCATTCGTG <sup>3</sup>

PAGE scheme

Schéma PAGE genotype

156/156 152/152 142/142 156/152 156/142 152/142







#### Used microsatellites

#### Examples:

*Pig* (10x) – S0068, S0107, SW24, SW840, SW353, SW936, SW353, S0070, SW72, TNFB

Cattle (11x) – TGLA53, TGLA122, TGLA126, TGLA227, ETH3, ETH10, ETH225, BM1824, BM2113, INRA23, SPS115

Horse (17x) – VHL20, AHT4, HTG4, HMS7, HTG6, AHT5, HMS6, ASB23, ASB2, HTG10, HTG7, HMS2, HMS3, ASB17, LEX3, HMS1, CA425

Dog (18x) – AHTk211, CXX279, INU055, REN169018, REN54P11, AHT137, AHTh260, AHTk253, INRA21, REN169D01, AHT121, AHTh171, FH2054, REN162C04, REN247M23, FH2848, INU005, INU30



raptor 5 MS

goat 14 MS

sheep 13 MS

• deer, cats, beavers, nutria, hens, geese, llamas




## Analysis procedure

- 1. Isolation of the genomic DNA of the tested individual
- 2. Multiplex PCR amplification of a panel of microsatellites
- 3. Fragmentation analysis using fluorescence capillary electrophoresis
- 4. Determination of peak sizes
- 5. Determination of alleles of individual microsatellites: discrimination by color and known range of alleles
- 6. Evaluation verification of parentage



# Multiplex PCR



Fluorescent colours(G5):

FAM - blue

VIC - green

NED – yellow

PET - red

LIZ – orange (size standard)

the allele size ranges of different microsatellites in the same color must not overlap




# Fragment analysis

Fluorescence capillary electrophoresis on a genetic analyzer - separation of microsatellites according to size






### Size standard





# Raw data - capillary electrophoresis








# Evaluation - genotypes

Example of one of the colours (6-FAM) and size standard (LIZ)





### **Partners:**





# Thank you for your attention!

This presentation has been supported by the Erasmus+ KA2 Cooperation Partnerships grant no. 2021-1-SK01-KA220-HED-000032068 "Innovation of the structure and content of study programs in the field of animal genetic and food resources management with the use of digitalisation - Inovácia obsahu a štruktúry študijných programov v oblasti manažmentu živočíšnych genetických a potravinových zdrojov s využitím digitalizácie". The European Commission support for the production of this presentation does not constitute an endorsement of the contents which reflects the views only of the authors, and the Commission cannot be held responsible for any use which may be made of the information contained therein.



Aleš Knoll



knoll@mendelu.cz



